

MTB Sampling Valve

- Lever-actuated on/off valve for storing and transporting samples in any GASTIGHT TLL syringe
- Kel-F hub needles attach to the valve
- Pressure tested to 100 psig; autoclavable
- Bore I.D., .039"



Ordering Information

Part Number	Description
86580	MTB Sampling Valve

Septum Adapter

- Store and sample syringe contents through the Septum Adapter
- Attaches to any TLL syringe or luer lock fitting
- Replaceable septum




Ordering Information

Part Number	Description
31335	Septum Adapter (TLL), includes 12 septa
75810	Replacement Septa, 7 mm dia., 12/pk

Removable (RN) Adapters

- For 250 µL and larger removable needle (RN) syringes
- Port I.D. of .039" allows sample spiking with a 701N syringe
- Provides connection to threaded fittings, tubing assemblies, tees and crosses



Ordering Information

Part Number	Description
35081	Female Luer/Removable Needle Hub (L)
35080	Male Luer/Removable Needle Hub (L)
35083	Male Luer Lock/Removable Needle Hub (L)